# NEW RELEASE

### Bringing Smiles to Everyday care

# The New Standard in Portable Toilets



Auto Wrapping Portable Toilet

## **o**mrabbou · **OKB**

Wrappon ORB - Seal away all your portable toilet concerns. Proudly Made in Japan.

## Simple to use!



1. Add the coagulant
It is used to solidify the moisture in waste.



#### 2. Press the button

After use, press the flush button to start the wrapping process. The process takes approximately 90 seconds to complete, and the remaining time is displayed on the remote control in a countdown format.



#### 3. Sealed wrap bag is ready

When a beep sound is heard, the wrapping process is complete. Pull out the tray and remove the wrapped bag after each use.

## Feature of wrappon ORB

#### Feature 1

#### Resin Body

The body is made of easy-to-clean resin.

#### Feature 2



#### Removable Armrests

The armrests can be easily removed to support smooth transfers.

#### Feature 3



#### Foldable Backrest

The backrest can be folded down easily using a rear lever.

#### Feature 4



#### **Bucket-Compatible**

An optional bucket is available as a backup for situations where power is not accessible.

#### Feature 5



#### Paper Holder

The toilet holder can be attached on either side.

#### Feature 6



#### Easy to Move

With the optional casters, the unit can be moved effortlessly.

#### Feature 7



#### Height-Adjustable Legs

The seat height can be adjusted to three levels: 40 cm, 43 cm, and 46 cm. Optional casters can also be attached for added convenience.

### MIHON SAFETY CO.,LTD.

#### http://www.wrappon.com

Wrappon Division, Nihon Safety Co., Ltd. 11F,Ichibancho Tokyu Bldg. 21 Ichibancho,Chiyoda-ku,Tokyo 102-0082 JAPAN E-mail:wrappon.usa@nihonsafety.com

**Qurappon** is a registered trademark of NIHONSAFETY CO.,LTD.

These product descriptions are effective as of February 2023. For product improvement purposes, the specifications and external appearance may be changed without prior notice.